





geovista

PRICE LIST 2014 - 2015



geovista



## NEW PRODUCTS

	Nucleotides	
GV-29001	dNTP Mixture each 2.5mM solution, 2.5mM each, Purity>99% (HPLC), 1ml	800
GV-29002	dNTP Mixture each 10mM solution , 10mM each, Purity>99% (HPLC), 1ml	3,250
GV-29003	dNTP Mixture each 25mM solution , 25mM each,Purity>99% (HPLC), 1ml.	6,450
GV-29004	dATP 100mMsolution, 100mM, Purity>99% (HPLC)	3,500
GV-29005	dCTP 100mMsolution, 100mM, Purity>99% (HPLC)	3,500
GV-29006	dGTP 100mM solution, 100mM, Purity>99% (HPLC)	3,500
GV-29007	dTTP 100mM solution, 100mM, Purity>99% (HPLC)	3,500
GV-29008	dUTP 100mM solution, 100mM, Purity>99% (HPLC)	3,500
	PCR/PCR Master Mix	
GV-29009	Taq DNA Polymerase, 5 U/ $\mu$ l, 1000U/pack is isolated from the E. coli cloning Thermus aquaticus.The molecular weight of the product is approximately 94 kDa. EasyTaq DNA polymerase has the 5to 3DNA polymerase activity and 5to 3 exonuclease activity without 3-5exonuclease activity.	2,000
GV-29010	Taq DNA Polymerase, 5 U/ $\mu$ l, 5000U/pack is isolated from the E. coli cloning Thermus aquaticus.The molecular weight of the product is approximately 94 kDa. EasyTaq DNA polymerase has the 5to 3DNA polymerase activity and 5to 3exonuclease activity without 3-5exonuclease activity.	8,000
GV-29011	Fast Taq DNA Polymerase, 2-9KB/min 5U/ $\mu$ l, 1000U/pack	5,375
GV-29012	Fast Taq DNA Polymerase, 2-9KB/min 5U/ $\mu$ l, 5000U/pack	21,500
GV-29013	Long Taq DNA Polymerase, upto 40kb, 5 U/ $\mu$ l, 1000U/pack	7,310
GV-29014	Hot Start DNA Polymerase, 5 U/ $\mu$ l, 1000U/pack	18,275
GV-29015	Pfu DNA Polymerase, 5 U/ $\mu$ l, 1000U/pack	11,825
GV-29016	M-MLV Reverse Transcriptase(RNase H-) 200U/ $\mu$ l, 1000U/pack	1,100
GV-29017	M-MLV Reverse Transcriptase(RNase H+) 200U/ $\mu$ l, 1000U/pack	1,100

GV-29018	This PCR master mix contains the Taq DNA Polymerase, dNTPs, and optimized buffer but without any dye in it. This Taq Mastermix is provided at 2× concentration and used at 1× concentration by adding template, primer, and H <sub>2</sub> O. It has the 5' to 3'DNA polymerase activity and 5' to 3' exonuclease activity without 3'-5' exonuclease activity. 200 rxn	16,340
GV-29019	SYBR Green Mix, without ROX, 2x conc. 400 rxn	15,000
GV-29020	SYBR Green Mix, with ROX, 2x conc. 400 rxn	15,000
GV-29021	Sybr Green I, 20X, for SYBR qPCR, 1ml	10,000
GV-29050	Powder Free Nytro Nitrile Powder-free Gloves medium 100/box, 8box/case	5,324
GV-29051	Powder Free Nytro Nitrile Powder-free Gloves large 100/box, 8box/case	5,324
GV-29052	Powder Free Nytro Nitrile Strong Powder-free Gloves small 100/box, 8box/case	7,000
GV-29053	Powder Free Nytro Nitrile Strong Powder-free Gloves medium 100/box, 8box/case	7,000
GV-29054	Powder Free Nytro Nitrile Strong Powder-free Gloves large 100/box, 8box/case	7,000
GV-29022	Test Tube Stands 12mm for 50 tubes, 2/pack	1,000
GV-29023	Multi Purpose Rack, for 0.5, 1.5, 2.0, 15, 50ml tubes, 4/case	1,500
GV-29100	Waste Bins, 2 liter, 5/case	1,750
GV-29101	Waste Bins, 5 liter, 2/case	1,200
GV-29102	Waste Bins, 15 liter, 1/case	1,200
GV-29024	Micro Centrifuge Racks, 1.5ml for 50 tubes, 5/case	1,900
GV-29025	Racks for PCR Tubes, 0.2ml for 96 tubes, 5/case	1,475
GV-29026	Cryo Cardboard Boxes, 1.8ml for 81 tubes, 10/case	1,300
GV-29027	Racks for PCR Tubes with lid, 0.2ml for 96 tubes, 5/case	950
GV-29028	Centrifuge Tube Rack, 15 & 50ml, 5/case	1,400

The background of the page features a complex, abstract molecular or crystal lattice structure composed of blue spheres and connecting lines, set against a dark blue gradient.

**Authorized Distributor**



**GEOVISTA**  
#18 & 19, 3rd Floor, Horamavu Agara Main Road  
Bangalore 560043- India, Tel: +91-80-32468095, Fax: +91-80-28443875  
E- Mail: [contact@geovista.co.in](mailto:contact@geovista.co.in), Web: [www.geovista.co.in](http://www.geovista.co.in)